


CryoCare Touch

Advanced sculpting precision



Intelligence. Confidence. Influence.



CryoCare Touch System

CryoCare Touch system is a fully integrated, highly precise ablation treatment system featuring a new intuitive touch screen interface, improved compact design, and software upgrades that offer a predictable therapy delivery with exceptional ease of technical operation.

Touch Screen Monitor

Includes an integrated touch screen user interface to control and monitor the procedure

Monitor Adjustment

May be tilted and rotated up to 190° for optimal viewing

Cart Handles

Handles can be locked in the up position for easy transport

Locking Wheels

2 lockable swiveling casters and 2 directional wheels

Real-time readout

Displays a real-time readout of the argon and helium system pressures during the procedure

Remote Control Keypad

Alternative control for operating the CryoCare Touch system

TempProbe Connections

Connect up to eight (A-H), independent TempProbes

Cryoprobe Connections

Connect up to eight (1-8), single-use, disposable cryoprobes

Compact Design

Compact, easy-to-operate console



Intelligence

Interventional radiologists can treat almost any cancer. At Varian, we're deeply committed to partnering with our customers in the development of innovative treatment solutions.

Confidence

The control, precision, and reliability offered by the CryoCare Touch system and V-Probe adjustable isotherm inspire more confident decision-making.

Influence

Varian is committed to offering products and services that facilitate collaboration and broaden the Interventional Radiologists' influential role in delivering enhanced patient outcomes.

Indications For Use

The CryoCare TOUCH system is intended for use in open, minimally invasive or endoscopic surgical procedures in the areas of general surgery, urology, gynecology, oncology, neurology, dermatology, ENT, proctology, pulmonary surgery and thoracic surgery. The system is designed to freeze/ablate tissue by the application of extreme cold temperatures including prostate and kidney tissue, liver metastases, tumors, skin lesions, and warts.

In addition, the CryoCare TOUCH system is intended for use in the following indications:

General Surgery

- Destruction of warts or lesions
- Palliation of tumors of the oral cavity, rectum and skin
- Ablation of breast fibroadenomas
- Ablation of leukoplakia of the mouth, angiomas, sebaceous hyperplasia, basal cell tumors of the eyelid or canthus area, ulcerated basal cell tumors, dermatofibromas, small hemangiomas, mucocele cysts, multiple warts, plantar warts, hemorrhoids, anal fissures, perianal condylomata, pilonidal cysts, actinic and seborrheic keratoses, cavernous hemangiomas, recurrent cancerous lesions

Urology

- Ablation of prostate tissue in cases of prostate cancer and benign prostatic hyperplasia

Gynecology

- Ablation of malignant neoplasia or benign dysplasia of the female genitalia

Oncology

- Ablation of cancerous or malignant tissue
- Ablation of benign tumors
- Palliative intervention

Neurology

- Freezing of nerve tissue in pain management/cryoanalgesia

Dermatology

- Ablation or freezing of skin

Proctology

- Ablation of benign or malignant growths of the anus or rectum
- Ablation of hemorrhoids

Thoracic Surgery

- Ablation of cancerous lesions

Achieve exceptional visualization and definitive control by pairing CryoCare Touch with V-Probe

Capable of producing an array of "ice-ball" sizes. Available in straight and right-angled configurations.

Exceptional Visualization^{1,2,3,4}

Clearly visualize the "ice-ball" on imaging in real time to confirm ablation zones adjacent to vulnerable structures (e.g., spinal column and the bowel) to minimize the risk of thermal injury. Overcome "cold sink" effects particularly centrally within the kidney.

Definitive Control³


Cytotoxic temperatures are delivered to the tissues within seconds due to rapid, high-pressure argon gas expansion by the Joule-Thomson effect. Varian cryoprobes are each equipped with an integrated temperature sensor to monitor probe performance during the procedure.

Adjustable Isotherms

The Varian V-Probe cryoprobe is capable of producing an array of "ice-ball" sizes. Control the size of the isotherm by moving an adjustable slider with your thumb. Available in straight and right angled configurations.


V-Probe Adjustable Cryoprobes

V-Probe cryoprobe offers 5 spherical "ice-ball" settings in a single probe. Available in both a right-angle and straight configuration.



V-Probe Isotherms

NOTE: The isotherms shown represent measurements collected using gelatin formula which approximates performance ($\pm 5\text{mm}$) in soft tissue at 100% gas for 10 minutes. The actual isotherm may vary.




Right Angle Cryoprobes

Right-angled design allows cryoprobe to be utilized in CT scanners. Available in multiple sizes and lengths.



Right Angle Cryoprobe Isotherms

NOTE: The isotherm measurements shown are for that particular probe model. Refer to the Directions for Use packaged with each probe for full isotherm measurement details. The isotherms shown represent measurements collected using gelatin formula which approximates performance ($\pm 5\text{mm}$) in soft tissue at 100% gas for 10 minutes. The actual isotherm may vary.



CryoCare System

Catalog Number	Product Description	Quantity
CRYO-38T	CryoCare Touch System: CryoCare Touch	1

V-Probe & V-Probe Kits

Catalog Number	Product Description	Shaft Diameter	Shaft Length	Quantity
CVA2400	2.4 mm V-Probe Cryoprobe	2.4 mm	15 cm	1
CRYO-207-V	Prostate Procedure Probe Pack – variable (6x CVA2400, 6x CRYO-55F, 1x CRYO-77)	2.4 mm	15 cm	1
CRYO-204-V	Prostate Procedure Probe Pack – variable (4x CVA2400, 4x CRYO-55F, 1x CRYO-77)	2.4 mm	15 cm	1

Right Angled V-Probe & Cryoprobes

Catalog Number	Product Description	Shaft Diameter	Shaft Length	Quantity
CVA2400RA	2.4 mm Right Angle V-Probe Cryoprobe	2.4 mm	15 cm	1
PCS-17R	1.7 mm Round Ice PerCryo Cryoprobe	1.7 mm	15 cm	1
PCS-17RS	1.7 mm Round Ice PerCryo Cryoprobe, Short	1.7 mm	7 cm	1
PCS-17	1.7 mm PerCryo Cryoprobe	1.7 mm	15 cm	1
RS-17	1.7 mm Round Ice PerCryo Cryoprobe	1.7 mm	15 cm	1
RS-17L	1.7 mm Long Shaft RenalCryo Cryoprobe	1.7 mm	23 cm	1
RS-24	2.4 mm RenalCryo Cryoprobe	2.4 mm	15 cm	1
RS-24L	2.4 mm Long Shaft RenalCryo Cryoprobe	2.4 mm	23 cm	1
PCS-24	2.4 mm PerCryo Cryoprobe	2.4 mm	15 cm	1
PCS-24L	2.4 mm Long Shaft PerCryo Cryoprobe	2.4 mm	23 cm	1
R3.8	3.8 mm RenalCryo CryoProbe	3.8 mm	13 cm	1
R3.8L	3.8 mm Long Shaft RenalCryo CryoProbe	3.8 mm	28 cm	1

Straight Cryoprobe

Catalog Number	Product Description	Shaft Diameter	Shaft Length	Quantity
CRYO-48-F	8.0 mm Cryoprobe	8 mm	29 cm	1

Temperature Probes & Other

Catalog Number	Product Description	Shaft Diameter	Shaft Length	Quantity
CRYO-55-F	TempProbe Device	16-gauge	14.5 cm	1

For a comprehensive list of available cryoprobes, please consult your Varian representative.

To order Varian cryoablation products, please contact Varian's Interventional Solutions Customer support team +1 877-467-2009 or vis.support@varian.com.

1. Zhong J, Wah TM. Renal ablation: current management strategies and controversies. Chin Clin Oncol. 2019;8(6):63.
2. Filippidis D, Mauri G, Marra P, Charalampopoulos G, Gennaro N, De cobelli F. Percutaneous ablation techniques for renal cell carcinoma: current status and future trends. Int J Hyperthermia. 2019;36(2):21-30.
3. Krokidis ME, Kitrou P, Spiliopoulos S, Karnabatidis D, Katsanos K. Image-guided minimally invasive treatment for small renal cell carcinoma. Insights Imaging. 2018;9(3):385-390.
4. Park SY, Won JY, Oh YT, Jung DC, Kim GM, Kim MD. Assessment of cold sink effect in postulated renal cryoablation by analyzing radiographic ice ball on computed tomography. Br J Radiol. 2018;20170951.

SAFETY STATEMENT: The CryoCare Touch System is intended for use in open, minimally invasive or endoscopic surgical procedures in the areas in general surgery, urology, gynecology, oncology, neurology (including cryoanalgesia), dermatology, ENT, proctology, pulmonary surgery and thoracic surgery. The system is designed to freeze/ablate tissue by the application of extreme cold temperatures, including prostate and kidney tissue, liver metastases, tumors, skin lesions, and warts. Use only as directed. As with all surgical procedures, the possibility of adverse reactions such as pain, fever, chills, sepsis, edema, perforation, ulceration, and hemorrhage may occur. For more information, please visit www.varian.com/safety.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Indications, contraindications, warnings, and instructions for use can be found in the product labeling supplied with each device.

varian

A Siemens Healthineers Company

varian.com/is

**USA, Corporate
Headquarters**

Varian Medical Systems
Palo Alto, CA 94304
Tel: 650.424.5700
800.544.4636

**USA, Varian's Interventional
Solutions Business Headquarters**

Varian Medical Systems
Austin, TX
Tel: 877.467.2009

**Legal
Manufacturer**

Varian Medical Systems
9825 Spectrum Drive
Building 2, Suite 275
Austin, TX 78717
USA

EC REP

**Authorized Representative
and Importer in the EU**

Varian Medical Systems
Nederland B.V.
Kokermolen 2
3994 DH Houten, The Netherlands
customer.relations@varian.com

vis.support@varian.com

CE
2797

Varian Medical Systems as a medical device manufacturer cannot and does not recommend specific treatment approaches. Specifications subject to change without notice. Not all features or products are available in all markets and are subject to change. Consult your Varian representative for country-specific product availability.

© 2023 Varian Medical Systems, Inc. VARIAN, CRYOCARE, PERCRYO, RENAL CRYO, TEMPPROBE and V-PROBE are trademarks of Varian Medical Systems, Inc., pending or registered U.S. Pat. & Tm. Off and other jurisdictions.