

DAKEWE

SurePrint C7 cassette printer

The compact allrounder



SurePrint C7 cassette printer

The future is non-contact laser marking

Enjoy smooth printing performance and savings on operational expenditure provided by non-contact UV laser printing technology

Compact and well-packed

Fits neatly into grossing station and labs

Continuous output

Special-engineered cassette output rack allows for non-stop printing, never stop the workflow

Intuitive handling and viewing of printed cassettes

Omnipotent to print on all sorts of cassettes

Capable to print all regular cassette configurations: Front or back hinged, with or without lid, one-piece or separate

UV laser tackles a wide range of polymers and prints well on different cassettes



Simple and powerful

Simple operation with friendly UI, yet powerful to adapt to various lab requirements

Easy integration with various LIS

Flexible template editing to provide layout fulfilling every need

Rich content including multiple barcode formats

Fast printing, long-lasting marks

Great barcode scannability thanks to 2,500 dpi high resolution

Permanent marks that resist processing solvents, heat and scratch

Hassle free maintenance

No more ink cartridges, ribbons, light bulbs or print heads

Eliminates the costs and downtime associating with these consumables

Care for the operators

Dedicated air purification system keeps the environment fresh

Smooth touch on all edges and corners



Specifications

Printing technology	Non-contact UV laser marking technology
Print content	Numbers, characters in multiple languages, symbols, 1D/2D barcodes, graphics
Print speed	Up to 2.5 sec/cassette
Cassette type	Cassettes with/without lid, one-piece cassettes; 35° - 45° angle
Air purification	Build-in air purification system with carbon filter
Load capacity	Up to 80 cassettes in 1 hopper
Output rack	Continuous output
System interface	LIS, HIS, PIC, PACS
Operating systems	Windows 7, Windows 10, Windows 11, Linux
Power	100 - 240 VAC, 50/60 Hz, 210 VA
Dimensions	320 × 380 × 440 mm
Weight	Approx. 25 kg



Sanova

Advancing health

Sanova Pharma GesmbH

Haidestraße 4, A-1110 Wien
Spitzermühlestraße 6, A-4713 Gallspach

Tel.: +43 1 801 04 28 01
Techn. Service: +43 1 801 04 25 87
medicalsystems@sanova.at
<https://medicalsystems.sanova.at>

