



GlideScope ClearFit

A Cover-Based Video Laryngoscope Solution



GlideScope[®] ClearFit[™]
verathon

GlideScope ClearFit Video Baton and Covers

GlideScope® ClearFit™ video laryngoscope addresses diverse airway management needs with versatile cover options, clear airway views, and compatibility across handheld and cart-based systems. It supports efficient and reliable airway management for both routine and difficult cases.



Versatility, Visualization, and Value in a Cover-Based Solution

Versatility for Clinician Preference

ClearFit helps you be prepared with a wide range of single-use covers. From traditional Mac and Miller to proven¹ Hyperangle covers, ClearFit accommodates diverse patient anatomies and clinician preference.

- 6 single-use cover options connect quickly with one reusable video baton
- Hyperangle cover design has a 22% lower profile compared to GVL® Stats²

Excellent Visualization for Clarity

Clearly view the patient's anatomy and quickly secure their airway. Elevate team-based care through visualization on large, high-resolution display options.

- Optimal camera placement for wide field of view
- Reveal™ anti-fog mechanism

Value for Cost Effective Care

Your fast-paced teams need efficient solutions to manage all types of airways. ClearFit, built with QuickConnect™ technology, enables interoperability across GlideScope handheld and cart-based systems, to streamline workflow and provide cost-effective patient care.

- Interchangeable use of 3.5" handheld and 10" or 15" cart-based monitors
- Reduced costs associated with inventory management, cleaning and training



1. Ruetzler K, Bustamante S, Schmidt MT, et al. Video Laryngoscopy vs Direct Laryngoscopy for Endotracheal Intubation in the Operating Room: A Cluster Randomized Clinical Trial. JAMA. 2024;331(15):1279–1286. doi:10.1001/jama.2024.0762.

2. Height at camera comparison of the ClearFit Hyperangle 4 cover to the GVL 4 Stat.

Designed for First-Pass Success

ClearFit™ video baton and covers deliver reliable GlideScope performance and flexibility for both routine and difficult airways. With a spring-loaded mechanism, you can seamlessly connect Mac, Miller or Hyperangle covers to suit your clinical preference and diverse patient anatomies, streamlining airway management and supporting first-pass success¹.



QuickConnect™ technology
facilitates seamless transition between
handheld and cart-based systems

Reusable baton design
built for repeated performance

Security tab
allows the single-use cover to
connect securely and release easily

Innovative retractable design
accommodates multiple cover styles

Single-use covers
minimize cross-contamination

Reveal anti-fog mechanism
resists lens fogging for
consistently clear airway views

6 cover styles
available for diverse
clinical preferences




Learn More
verathon.com/clearfit

1. Ruetzler K, Bustamante S, Schmidt MT, et al. Video Laryngoscopy vs Direct Laryngoscopy for Endotracheal Intubation in the Operating Room: A Cluster Randomized Clinical Trial. JAMA. 2024;331(15):1279–1286. doi:10.1001/jama.2024.0762.

GlideScope Portfolio Compatibility for Standardization

ClearFit™ video laryngoscope offers flexibility for routine and difficult airway management. Its innovative design allows your team to choose the monitor size and the cover style that best fits specific clinical needs or preferences. Additionally, the magnetic QuickConnect™ technology drives workflow efficiency and integrates the GlideScope® ecosystem—offering compatibility with a broad range of handheld and cart-based monitors, ClearFit, and Spectrum™ QC single-use video laryngoscopes, which include pediatric size and style options.

Compatible Monitors



GlideScope® Core™ 15



GlideScope® Core™ 10




GlideScope® Core™
2 m QuickConnect™ Cable



GlideScope® Go™ 2
Portable, 3.5".

Works with GlideScope Core 2m QuickConnect Cable and GlideScope Core workstation.

GlideScope® ClearFit™ Video Baton with Single-Use ClearFit Covers




Hyperangle 3, 4

Miller 2

MAC 2, 3, 4

GlideScope® Spectrum QC eco™ Single-Use Video Laryngoscopes



Hyperangle S1, S2, S2.5, S3, S4

Miller S00, S0, S1, S2

MAC S1, S2, S3, S4



See the difference up close.
Request a demo today.

Single-use covers and blades are sold in boxes of 10.

verathon

verathon.com

Verathon Inc.
20001 North Creek Parkway
Bothell, WA 98011
USA
Tel: +1 800 331 2313
(USA & Canada only)
Tel: +1 425 867 1348
Fax: +1 425 883 2896

EC REP

Verathon Medical (Europe) B.V.
Willem Fenengastraat 13
1096 BL Amsterdam
The Netherlands
Tel: +31 (0) 20 210 30 91
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited
Unit 9, 39 Herbert Street
St Leonards NSW 2065
Australia
Tel: 1800 613 603 (Within Australia)
Fax: 1800 657 970
Tel: +61 2 9431 2000 (International)
Fax: +61 2 9475 1201