

PRODUCT INFO

INSTRUMENTS, CHEMICALS AND CONSUMABLES FOR CELL DIAGNOSTICS

PRODUCT GROUP CHEMICALS

HISTOLAB-CLEAR

Histolab-Clear is a solvent, consisting of ISO-paraffins and replaces Xylen both at the dehydration and mounting process. The advantage of Histolab-Clear is that it has a higher flash point and boiling point than Xylen.

Histolab-Clear is an almost odourless, non-carcinogenic solvent which have the properties to solve paraffin. Histolab-Clear minimizes the risk of getting to hard tissues at dehydration, particularly tissues containing keratin.

Each package is marked with batch number for traceability.



5 L
ART CODE: 14250

SPECIFICATIONS SOLVENT HISTOLAB-CLEAR

Art code:	14250	5 L
	14275	200 L
	14277	400 L in container

Shelf life: Minimum 5 years



Histolab Products AB
Södra Långebergsgatan 36
SE-421 32 Gothenburg, Sweden
Tel: +46 31 709 30 30
Fax: +46 31 709 30 40
mail@histolab.se
www.histolab.se