



## SteriStay Extension – extend your protective instrument table

The SteriStay Extension increases your sterile area even further than the regular instrument table. The SteriStay Extension is an attachment that easily attaches on to the existing SteriStay instrument table. Once not needed, it is just as easy to remove it.

SteriStay Extension ensures that instruments and implants remain sterile during both preparation and the entire surgical procedure. It guarantees an ultraclean air environment with  $<5$  cfu/m<sup>3</sup> air around your sterile instruments at all times.



*“Everything which is to come into contact with the wound has been made sterile, except the air, which is in contact with everything” Whyte 1973*

SteriStay Extension is used in the same way as SteriStay instrument table. It can be used in all types of operating rooms and surgical preparation rooms, independent of ventilation system, where sterile instruments are handled. By using SteriStay Extension you can expand the ultra clean zone even further in the operating room. Setting up the table in advance under safe and secure hygienic conditions minimizes setup time between surgeries.

SteriStay Extension is an accessory to SteriStay that bears the CE marking and is a Class I device according to the MDD (Medical Device Directive) 93/42 EEC.



**Intertek**

Toul Meditech's Quality Management System has been assessed and registered by Intertek Certification AB as conforming to the requirements of ISO 13485:2012.

SteriStay is protected by patent.

## Specifications

| <b>Physical Data:</b> |   |
|-----------------------|---|
| Size and weight:      | Size: L 40 x W 60<br>Total work surface: 173 cm<br>Weight: 5 kg |
| Load limitation:      | Extension may be loaded with max 5 kg                           |
| <b>Material:</b>      |   |
| Surface:              | ABS plastic   |
| Fitting frame:        | Stainless steel   |
| <b>Attachment:</b>    |   |
| Assemble              | Snap lock   |

Toul Meditech AB  
Tunbytorpsgatan 31  
SE 721 37 Västerås, Sweden

Phone: int + 46 21 13 50 00  
Fax: int + 46 21 13 86 45

Internet: [www.toulmeditech.com](http://www.toulmeditech.com)  
E-mail: [info@toulmeditech.com](mailto:info@toulmeditech.com)