

# UPLIFT

Das präzise **laparoskopische Sakrokolpopexie** System

## SPEZIELLES DESIGN FÜR DIE PERFEKTE JUSTIERUNG

Einfache Technik

Einfache Ankerfixierung

Laparoskopischer Ankerapplikator der neuesten Generation


35 cm Schaft  
Ø 5mm.

Unabhängige anteriore und  
posteriore Spannungsregelung

## EINFACHES UND SCHNELLES SYSTEM FÜR DIE LAPAROSKOPISCHE SACROCOLPOPEXY

Makroporöses, weiches und vorgeschchnittenes Netz aus 100 % monofilen Polypropylen


- ZEITERSPARNIS GEGENÜBER LAPAROSKOPISCHEM NÄHEN
- PRÄZISE KONTROLLE DER ANKERPLATZIERUNG
- ATRAUMATISCHE ANKER



## GLEICHES MESH MIT 2 UNTERSCHIEDLICHEN PORENGRÖßEN



**NEOMESH**  
UPLIFT  
KITUPL-05



**NEOPORE**  
UPLIFT04  
KITUPL04-05

ARTIKELNUMMER	PRODUKT	BESCHREIBUNG
UPLIFT	UPLIFT	Neomesh Uplift Mesh für laparoskopische Sakrokolpopexie
KITUPL-05	Kit UPLIFT	Neomesh Uplift Mesh für laparoskopische Sakrokolpopexie Inhalt: 1 Neomesh Uplift Mesh, 1 Uplift Applikator + 5 Uplift Anker
UPLIFT04	UPLIFT	Neopore Uplift Mesh für laparoskopische Sakrokolpopexie
KITUPL04-05	Kit UPLIFT	Neopore Uplift Mesh für laparoskopische Sakrokolpopexie Inhalt: 1 Neopore Uplift Mesh, 1 Uplift Applikator + 5 Uplift Anker
UPLANCHOR-05	UPLIFT ANKER	5 Uplift Anker
UPLAPLIC 5A	UPLIFT ANKER SYSTEM	1 Uplift Applikator + 5 Uplift Anker

info@neomedic.com



The Continence Company